



Assay	What it Detects	What it Doesn't Detect	TAT ⁱ	When to Use
G-Banded Karyotyping	<ul style="list-style-type: none"> Microscopic genomic abnormalities (>5-10 Mb) <ul style="list-style-type: none"> Inversions Duplications/deletions Balanced and unbalanced translocations Aneuploidies >10% mosaicism (for example: cultures where >1 of 10 cells are <u>trisomy 12</u>) 	<ul style="list-style-type: none"> Submicroscopic genomic abnormalities (<5 Mb) 	7-10 days	<ul style="list-style-type: none"> As a baseline genomic screen <ul style="list-style-type: none"> At derivation of cell lines At the start of experimental protocols To assess and monitor genomic stability (for example: every 5-10 passages of cell culture) At conclusion of experiments (prior to publication) For cell line banking When publication-quality karyotypes are needed
Spectral Karyotyping (SKY)	<ul style="list-style-type: none"> Microscopic genomic abnormalities (>5-10 Mb) <ul style="list-style-type: none"> Balanced and unbalanced translocations Aneuploidies 	<ul style="list-style-type: none"> Submicroscopic abnormalities (<5 Mb) Inversions Duplications/deletions 	10-15 days	<ul style="list-style-type: none"> As an adjunct to g-banded karyotyping <ul style="list-style-type: none"> To define complex rearrangements To identify marker chromosomes When publication-quality spectral karyotypes are needed
Microarray by SNP and/or CGH	<ul style="list-style-type: none"> Microscopic genomic abnormalities (>5-10 Mb) Genomic gains and losses (>50 Kb) <ul style="list-style-type: none"> Copy number variants Duplications/deletions Unbalanced translocations Aneuploidies >~10% mosaicism (for example: cultures where >1 of 10 cells are trisomy 12) Loss of Heterozygosity (LOH) / Absence of Heterozygosity (AOH) - detected by SNP 	<ul style="list-style-type: none"> Balanced translocations Inversions <~10% culture mosaicism (for example: cultures where 1 of 10 cells is trisomy 12) Chromosomal position of genomic gains Loss of Heterozygosity (LOH) / Absence of Heterozygosity (AOH) - not detected by CGH alone 	14-21 days	<ul style="list-style-type: none"> As a baseline genomic screen <ul style="list-style-type: none"> To detect submicroscopic (<5 Mb) abnormalities To identify amplified or deleted genes of interest In conjunction with g-banded karyotyping <ul style="list-style-type: none"> To define unbalanced translocation breakpoints For research of genomic copy number change <ul style="list-style-type: none"> To identify structural variation within populations or disease cohorts To develop a cell line copy number variant profile
Fluorescence In Situ Hybridization (FISH)	<ul style="list-style-type: none"> Genomic sequence of interest <ul style="list-style-type: none"> Duplications or deletions >20 Kb >2% mosaicism (for example: cultures where >2 of 100 cells are trisomy 12) Chromosomal location of genomic gains Chromosome fusions (breakaparts) 	<ul style="list-style-type: none"> Changes in regions other than the probe-specific sequence 	10-15 days	<ul style="list-style-type: none"> To confirm findings and refine breakpoints detected by g-banded karyotyping To confirm findings and localize genomic gains detected by aCGH analysis As a screen for microdeletions/duplications of known targets
Fast FISH	<ul style="list-style-type: none"> Greater than 50% of recurrent trisomies in both human ES and iPS cell lines.ⁱⁱ <ul style="list-style-type: none"> Trisomies of human chromosomes X, 8, 12 and 17. 	<ul style="list-style-type: none"> Aneuploidies of chromosomes other than X, 8, 12 and 17 Structural chromosome abnormalities 	1-2 days	<ul style="list-style-type: none"> As a rapid screen for large numbers of clones As a cost-effective tool for monitoring aneuploidy
Short Tandem Repeat Analysis (STR)	<ul style="list-style-type: none"> STR polymorphisms for 15 loci plus amelogenin (Promega® PowerPlex® 16) Probability of matching identity to an existing STR profile 	<ul style="list-style-type: none"> STR polymorphisms in areas other than those represented in Promega® PowerPlex® 16 	10-15 days	<ul style="list-style-type: none"> To monitor identity of a cell line To confirm relationship of iPS cells to their parent line To establish an STR profile of a newly derived or reprogrammed cell line To rule out culture cross-contamination

ⁱ Turn-around times provided are based on provision of sufficient mitotically active hES or iPS cultures grown in matrigel/TeSR or MEF/hES media conditions.

ⁱⁱ Based on n=1220 hES/iPS cell cultures karyotyped by WiCell Cytogenetics Laboratory.