

Thaw and Culture Details

Cell Line Name	MIN31i-33363.D.3C2			
WiCell Lot Number	WB67625			
Provider	Massachusetts General Hospital			
Banked By	WiCell			
Thaw and Culture Recommendations	WiCell recommends thawing 1 vial into 3 wells of a 6 well plate using mTeSR™Plus and Matrigel®.			
Protocol	WiCell Feeder Independent Pluripotent Stem Cell Protocol			
Culture Platform Prior to Freeze	Feeder Independent			
	Medium: mTeSR™Plus			
	Matrix: Matrigel®			
Passage Number	p24 These cells were cultured for 23 passages prior to freeze and post reprogramming. WiCell adds +1 to the passage number at freeze to best represent the overall passage number of the cells at thaw. Plated cells at thaw should be labeled passage 24.			
Date Vialed	29-January-2021			
Vial Label	MIN31i-33363.D.3C2 p24 WB67625			
Biosafety and Use Information	Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans.			

Testing Performed by WiCell

Test Description	Test Provider	Test Method	Test Specification	Result
Karyotype by G-banding	WiCell	SOP-49	Expected karyotype	See Report
Post-Thaw Viable Cell Recovery	WiCell	SOP-99	≥ 15 Undifferentiated Colonies prior to passage, ≤ 30% Differentiation prior to passage, and recoverable attachment after passage	Pass
Identity by STR	UW Translational Research Initiatives in Pathology Laboratory	PowerPlex 16 HS System by Promega	Defines STR profile of deposited cell line	Pass
Sterility	Steris	ST/07	Negative	Pass
Mycoplasma	WiCell	SOP-79	Negative	Pass

Approval Date	Quality Assurance Approval
10-March-2021	3/10/2021 X JKG IKG Quality Assurance Signed by Gay, Jenna



Chromosome Analysis Report: 084804

Date Reported: Monday, February 22, 2021


Cell Line: MIN31i-33363.D.3C2-WB67625

Submitted Passage #: 24

Date of Sample: 2/11/2021

Specimen: Human IPSC

Results: 46,XY



Cell Line Sex: Male

Reason for Testing: none given

Investigator: WiCell Stem Cell Bank, WiCell

Cell: 35

Slide: G02

Slide Type: Karyotype

Total Counted: 20
Total Analyzed: 8

Total Karyogrammed: 4

Band Resolution: 425 - 475

QC Review By: _

Interpretation:

Date:

This is a normal karyotype; no clonal abnormalities were detected at the stated band level of resolution.

Sent By:____ Sent To:_

cell populations in this specimen (i.e.,mosaicism) is limited by the number of metaphase cells examined, documented here as "# of cells counted".

Completed by:	, CG(ASCP)
Reviewed and Interpreted by:	, PhD, FACMG

Limitations: This assay allows for microscopic visualization of numerical and structural chromosome abnormalities. The size of structural abnormality that can be detected is >3-10Mb, dependent upon the G-band resolution obtained from this specimen. For the purposes of this report, band level is defined as the number of G-bands per haploid genome. It is documented here as "band level", i.e., the range of bands determined from the four karyograms in this assay. Detection of heterogeneity of clonal

This assay was conducted solely for listed investigator/institution. The results of this assay are for research use only. Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.



Short Tandem Repeat

Requestor: WiCell Stem Cell Bank, WiCell Samples Received: 11Feb21 STR Amplification Date: 15Feb21

	MIN31i-
Sample Name	33363.D.3C2-
	WB67625 p.24
Label on tube	84804
FGA	
TPOX	
D8S1179	
vWA	
Amelogenin	Identifying
Penta_D	information has
CSF1PO	been redacted to
D16S539	protect donor confidentiality. If
D7S820	more information
D13S317	is required,
D5S818	please contact info@wicell.org
Penta_E	ii ii o ee wiccii.org
D18S51	
D21S11	
TH01	
D3S1358	
Allelic Polymorphisms	28
Matches*	
Comments	

*Note: The STR profile of the following sample is an exact match for the given sample/samples.



Short Tandem Repeat

Form SOP-89.01 Version 3.0

Requestor: WiCell Stem Cell Bank, WiCell Samples Received: 11Feb21 STR Amplification Date: 15Feb21

Results: The genotypic profiles comprise a range of 28 allelic polymorphisms across the 15 STR loci analyzed.

Interpretation: The concentration of DNA required to achieve an acceptable STR genotype (signal/noise) was equivalent to that required for the standard procedure (~1 ng/amplification reaction) from human genomic DNA. These results suggests that the cells submitted correspond to the cell lines as named and were not contaminated with any other human cells or a significant amount of mouse feeder layer cells.

Sensitivity: Sensitivity limits for detection of STR polymorphisms unique to either this or other human cell lines is ~2-5%.

	2/17/2021	2/17/2021	2/18/2021
X		X	X
Tech #1 Characterization Signed by:		Tech #2 Characterization Signed by:	OA Review Quality Assurance Signed by:

Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.

Raw data is available upon request.

Native Product Sterility Report



SAMPLE #:

21020702

DATE RECEIVED:

11-Feb-21

TEST INITIATED:

12-Feb-21

TEST COMPLETED:

26-Feb-21

SAMPLE NAME / DESCRIPTION:

504 S Rosa Road, Rm 101

Madison, WI 53719

MIN28i-35833.A-WB67616

MIN31i-33363.D.3C2-WB67625

MIN29i-35833.B-WB67612

MIN30i-33109.2G-WB67613 MIN27i-35326.K-WB67617

JHU042i-WB67618

MIN33i-33109.2G.1A4-WB67621

WA09-WB67619

MIN32i-33109.2B.3A12-WB67622 MIN34i-33109.2G.2F2-WB67623

UNIQUE IDENTIFIER:

N/A

TEST RESULTS:

WiCell

'roudh' Oantin	
3	Growth) - Contro O 2 Negative

TEST SUMMARY:

# Samples	Media Type	Volume (mL)	Incubation Temperature (° C)	Incubation Duration (Days)
10	TSB	40	20-25	14
10	FTG	40	30-35	44

REFERENCE:

Processed according to LAB-003: Sterility Test Procedure

PD #:

000053

TEST METHODOLOGY:

ISO - Direct Transfer

COMMENTS:

NA

REVIEWED BY

DATE OF MAR 2021

Specific test results may not be indicative of the characteristics of any other samples from the same lot or similar lots. This test report shall not be reproduced, except in full, without prior written approval. Liability is limited to the costs of the tests. Results applied to samples as received.



Mycoplasma Assay Report

PCR-based assay performed by WiCell
WiCell
16Feb21

Sample Name	Result	Interpretation
PENN005i-35-3-DB36317 p.14 (84799)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
PENN006i-149-1-DB36519 p.13 (84800)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
MIN31i-33363.D.3C2-WB67625 p.24 (84804)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
INC149 11FEB21AP (84805)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
INC169 11FEB21MMM (84806)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
INC123 11FEB21KR 1 OF 2 (84808)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
INC123 11FEB21KR 2 OF 2 (84836)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
Positive (+) Control	Positive	
Negative (-) Control	Negative	

Reported by: , Assistant Research Specialist Reviewed by: , Assistant Research Specialist

Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.

A gel image is available upon request.