

### **Thaw and Culture Details**

Cell Line Name	MCW110i-U2170			
WiCell Lot Number	WB67204			
Parent Material	MCW110i-U2170-DB66418			
Provider	Medical College of Wisconsin – Laboratory of Dr. Ulrich Broeckel			
Banked By	WiCell			
Thaw and Culture Recommendations	WiCell recommends thawing 1 vial into 3 wells of a 6 well plate.			
Culture Platform	Feeder Independent			
	Medium: TeSR <sup>™</sup> -E8 <sup>™</sup>			
	Matrix: Matrigel®			
Protocol	WiCell Feeder Independent E8 Medium Protocol			
Passage Number	p16 These cells were cultured for 15 passages prior to freeze and post colony selection. WiCell adds +1 to the passage number at freeze to best represent the overall passage number of the cells at thaw. Plated cells at thaw should be labeled passage 16.			
Date Vialed	25-May-2019			
Vial Label	MCW110i-U2170 p16 WB67204			
Biosafety and Use Information	Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans.			

## **Testing Performed by WiCell**

Test Description	Test Provider	Test Method	Test Specification	Result	
	WiCell	SOP-CH-003	Expected karyotype	See Report	
Karyotype by G-banding	<b>Results:</b> 47,XY, +12[2]/46,XY[18 <i>Interpretation:</i> This is an abnormal karyotype. An extra copy of chromosome 12 (trisomy 12 present in two of twenty cells examined. Gain of chromosome 12 is recurrently acquired in pluripotent stem cell cultures. No other clonal abnormalities were detected at the stated band level of resolution.				
Post-Thaw Viable Cell Recovery	WiCell	SOP-CH-305	<ul> <li>≥ 15 Undifferentiated Colonies prior to passage,</li> <li>≤ 30% Differentiation prior to passage, and recoverable attachment after passage</li> </ul>	Pass	
Identity by STR	UW Translational Research Initiatives in Pathology Laboratory	PowerPlex 16 HS System by Promega	Defines STR profile of deposited cell line	Pass	
Sterility	Steris	ST/07	Negative	Pass	
Mycoplasma	WiCell	SOP-CH-044	Negative	Pass	

©2020 WiCell Research Institute

The material provided under this certificate has been subjected to the tests specified and the results and data described herein are accurate based on WiCell's reasonable knowledge and belief. Appropriate Biosafety Level practices and universal precautions should always be used with this material. For clarity, the foregoing is governed solely by WiCell's Terms and Conditions of Service, which can be found at http://www.wicell.org/privacyandterms.



# **Testing Reported by Provider**

The Provider stated that some or all of the additional analyses listed below may have been performed for this cell line. For more information, publication and dbGaP links, where available, are provided on the cell line specific web page on the WiCell website.

- Tra1-60 marker expression
- mRNA expression by qPCR
- Infinium<sup>®</sup> Expanded Multi-Ethnic Genotyping Array (MEGA<sup>EX</sup>)

Approval Date	Quality Assurance Approval
18-June-2020	6/18/2020 XG Quality Assurance Signed by Gay, Jenna

©2020 WiCell Research Institute




Date Reported: Monday, November 25, 2019 Cell Line: MCW110i-U2170-WB67204 15153 Passage#: 16 Date of Sample: 11/20/2019 Specimen: Human IPSC Results: 47,XY,+12[2]/46,XY[18]

Cell Line Sex: Male Reason for Testing: Lot release

WiCell

Investigator:



Cell: 26 Slide: G04 Slide Type: Karyotype Total Counted: 20 Total Analyzed: 8 Total Karyogrammed: 4 Band Resolution: 450 - 475

#### Interpretation:

This is an abnormal karyotype. An extra copy of chromosome 12 (trisomy 12) is present in two of twenty cells examined. Gain of chromosome 12 is recurrently acquired in pluripotent stem cell cultures.

No other clonal abnormalities were detected at the stated band level of resolution.

Completed by: Reviewed and Interpreted by:		, Ph.D.		
Date:	Sent By:	Sent To:	QC Review By:	

Limitations: This assay allows for microscopic visualization of numerical and structural chromosome abnormalities. The size of structural abnormality that can be detected is >3-10Mb, dependent upon the G-band resolution obtained from this specimen. For the purposes of this report, band level is defined as the number of G-bands per haploid genome. It is documented here as "band level", i.e., the range of bands determined from the four karyograms in this assay. Detection of heterogeneity of clonal cell populations in this specimen (i.e., mosaicism) is limited by the number of metaphase cells examined, documented here as "# of cells counted".

This assay was conducted solely for listed investigator/institution. The results of this assay are for research use only. Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.



# Short Tandem Repeat Analysis



characterization@wicell.org (608) 316-4145

Receive Date: 12/05/19 Report Sent: 12/16/19 Assay Date: 12/10/19 File Name: STR 191212 wmr Report Date: 12/16/19

Department of Pathology and Laboratory Medicine TRIP Laboratory (Molecular) https://research.pathology.wisc.edu/trip-home/ (608) 265-9168

#### Sample Report:

15153-STR Sample Name on Tube: 15153-STR 139.4 ng/μL, (A260/280=1.90) Sample Type: Cells Cell Count: ~2 million cells Requestor: WiCell Research Institute Quality Assurance Department

STR Locus	STR Genotype Repeat #	STR Genotype
FGA	16-18,18.2,19,19.2,20,20.2,21,21.2,22, 22.2, 23, 23.2, 24, 24.2, 25, 25.2, 26-30, 31.2, 43.2,	Identifying
TROM	<u>44.2,45.2, 46.2</u> <u>6-13</u>	information has
TPOX	ND37-12	been redacted to
D8S1179	7-18	protect donor
vWA	10-22	confidentiality. If
Amelogenin	X,Y	more information
Penta D	2.2, 3.2, 5, 7-17	is required,
CSF1PO	6-15	please, contact
D16S539	5, 8-15	- <u>WiCell's Technical</u> Support.
D7S820	6-14	
D13S317	7-15	
D5S818	7-16	
Penta E	5-24	
D18S51	8-10, 10.2, 11-13, 13.2, 14-27	
D21S11	24,24.2,25,25.2,26-28,28.2,29,29.2, 30, 30.2,31, 31.2,32,32.2,33,33.2, 34,34.2,35,35.2,36-38	
<b>TH01</b>	4-9,9.3,10-11,13.3	
D3S1358	12-20	

<u>Results:</u> Based on the 15153-STR cells submitted by WiCell QA dated and received on 12/05/19, this sample (Label on Tube: 15153-STR) defines the STR profile of the human cell line MCW110i-U2170 comprising 25 allelic polymorphisms across the 15 STR loci analyzed.

<u>Interpretation:</u> No STR polymorphisms other than those corresponding to the human MCW110i-U2170 cell line were detected and the concentration of DNA required to achieve an acceptable STR genotype (signal/ noise) was equivalent to that required for the standard procedure (~1 ng/amplification reaction) from human genomic DNA. This result suggests that the 15153-STR sample submitted corresponds to the MCW110i-U2170 cell line and was not contaminated with any other human cells or a significant amount of mouse feeder layer cells.

<u>Sensitivity:</u> Sensitivity limits for detection of STR polymorphisms unique to either this or other human cell lines is ~2-5%.

X RMB	Digitally Signed on 12	:/16/19	X WMR	Digitally Signed on	12/16/19
TRIP La	, BA boratory, Molecular		UWHC Mole	, PhD, Director / Co-Direct cular Diagnostics Laboratory / UW	

Testing was accomplished by analysis of human genetic polymorphisms at STR loci. This methodology has not vet been approved by the FDA and is for investigational use only. Acknowledge TRIP in your publications, posters & presentations. For details, see: https://research.pathology.wisc.edu/acknowledging-trip/ Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at https://www.wicell.org/media.acux/ca76d97c-862a-43f3-b02a-ab2d1e619100. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.

# Native Product Sterility Report



	SAMPLE #:	19060913
WiCell	DATE RECEIVED:	12-Jun-19
504 S Rosa Road, Rm 101	TEST INITIATED:	14-Jun-19
Madison, WI 53719	TEST COMPLETED:	28-Jun-19

SAMPLE NAME / DESCRIPTION:	STAN204i-448C1 WB67189 14791
	MCW013i-A2767 WB67191 14792
	JHU242i DB37058 14793
	MCW085i-40002118 WB67193 14794
	MCW081i-U7128 WB67194 14795
	STAN043i-124-1 WB67196 14796
	STAN038i-118-2 WB67197 14797
	MCW007i-U2456 WB67198 14798
	MCW096i-40000169 WB67199 14799
	MCW074i-40002460 WB67203 14800
	MCW110i-U2170 WB67204 14801
	STAN044i-124-2 WB67206 14802
	MCW105i-U2130 WB67207 14803
	MCW103i-40000237 WB67208 14804
	MCW101i-40001005 WB67209 14805
	hIPSC-Di21-c2-4-4 WB67210 14806
	WA07 WB67212 14807
	WA07 WB67213 14808
	MCW021i-50001743 WB67214 14809
	hIPSC-Di21-c2-4-3 WB67215 14810

#### **UNIQUE IDENTIFIER:**

TEST RESULTS:	# Tested	# Positives (Growth)	- Control
	20	0	2 Negatives

NA

TEST	SU	MMA	RY:
------	----	-----	-----

TEST SUMMARY:	# Samples	Media Type	Volume (mL)	Incubation Temperature (° C)	Incubation Duration (Days)
	20	TSB	40	20-25	14
	20	FTG	40	30-35	14
REFERENCE:		Processed accord	ding to LAB-003: St	erility Test Procedu	ire

PD #:

TEST METHODOLOGY:

000053 USP - Direct Transfer

# Native Product Sterility Report



COMMENTS:

NA

REVIEWED BY

DATE 28JUNIA

Specific test results may not be indicative of the characteristics of any other samples from the same lot or similar lots. This test report shall not be reproduced, except in full, without prior written approval. Liability is limited to the costs of the tests. Results applied to samples as received.

N



### Mycoplasma Assay Report

PCR-based assay performed by WiCell WiCell 25Nov19

Sample Name	Result	Comments/Suggestions
MCW110i-U2170-WB67204 15153 (79186)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
MCW107i-40000886-WB67227 15154 (79187)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
MCW101i-40001005-WB67209 15159 (79188)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
WC060i-226-1-2-22-WB67334 15142 (79189)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
Positive (+) Control	Positive	
Negative (-) Control	Negative	

Reported by: Cell Culture Specialist Reviewed by: Cell Culture Specialist

Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.

A gel image is available upon request.