



Thaw and Culture Details

Cell Line Name	CREM043i-SA210-1
WiCell Lot Number	DB51680
Provider	Boston University – Laboratory of Dr. Martin Steinberg
Banked By	Boston University - Laboratory of Dr. Gustavo Mostoslavsky
Thaw and Culture Recommendations	The Provider recommends thawing 1 vial into 2 wells of a 6 well plate.
Culture Platform	Feeder Independent
	Medium: mTeSR™1
	Matrix: Matrigel®
Protocol	WiCell Feeder Independent mTeSR™1 Protocol
Passage Number	p6 These cells were cultured for 6 passages prior to freeze and post colony picking. Therefore, plated cells at thaw should be labeled passage 7.
Date Viald	02-July-2013
Vial Label	hiPSC SA210-1 P6 7.2.13 AS
Biosafety and Use Information	Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans.

Testing Reported by Provider

The Provider stated that some or all of the additional analyses listed below may have been performed for this cell line. For more information, publication and dbGaP links, where available, are provided on the cell line specific web page on the WiCell website.

- Digital Genome Sequencing
- Infinium® Expanded Multi-Ethnic Genotyping Array (MEGA^{EX})

Please note: Prior to shipment of these cells, WiCell will perform the following characterization assays: post-thaw viable recovery, identity by STR, sterility, mycoplasma, and karyotype.

Approval Date	Quality Assurance Approval
25-June-2018	6/25/2018  RK RK Quality Assurance Signed by: Kremers, Erik