



Thaw and Culture Details

Cell Line Name	CREM005i-SS2-1GAG
WiCell Lot Number	DB66769
Provider	Boston University – Laboratory of Dr. Martin Steinberg
Banked By	Boston University - Laboratory of Dr. Gustavo Mostoslavsky
Thaw and Culture Recommendations	The Provider recommends thawing 1 vial into 2 wells of a 6 well plate.
Culture Platform	Feeder Independent
	Medium: mTeSR™1
	Matrix: Matrigel®
Protocol	WiCell Feeder Independent mTeSR™1 Protocol
Passage Number	p39 These cells were cultured for 39 passages prior to freeze and post colony picking. Therefore, plated cells at thaw should be labeled passage 40.
Date Viald	19-December-2016
Vial Label	SS2-1 GAG p39 hiPSC/mTeSR 12/19/16 SMP
Biosafety and Use Information	Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans.

Testing Reported by Provider

The Provider stated that some or all of the additional analyses listed below may have been performed for this cell line. For more information, publication and dbGaP links, where available, are provided on the cell line specific web page on the WiCell website.

- Digital Genome Sequencing
- Infinium® Expanded Multi-Ethnic Genotyping Array (MEGA^{EX})

Please note: Prior to shipment of these cells, WiCell will perform the following characterization assays: post-thaw viable recovery, identity by STR, sterility, mycoplasma, and karyotype.

Approval Date	Quality Assurance Approval
25-June-2018	<p style="text-align: right;">6/25/2018</p> <p>X RK</p> <p><small>RK Quality Assurance Signed by Kremers, Erik</small></p>