



Certificate of Analysis

Thaw and Culture Details

Cell Line Name	SCR4301i	
WiCell Lot Number	DB42912	
Provider/Client	The Scripps Research Institute – Laboratory of Dr. Eric Topol	
Banked By	Gladstone Institutes – Laboratory of Dr. Sheng Ding	
Thaw and Culture Recommendations	WiCell recommends thawing 1 vial into 4 wells of a 6 well plate using mTeSR™ 1 and Matrigel®.	
Protocol	WiCell Feeder Independent Pluripotent Stem Cell Protocol	
Culture Platform Prior to Freeze	Medium: mTeSR™ 1	Matrix: Matrigel®
Passage Number	p10 Cells were cultured for 10 passages prior to freeze and post colony selection . Plated cells at thaw should be labeled passage 11.	
Date Vial	18-MAY-2016	
Vial Label	C00793, Passage 10, May-18-2016	
Biosafety and Use Information	Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans.	



Certificate of Analysis

Results

Test Description	Test Provider	Test Method	Test Specification	Result
Karyotype	WiCell	G-T-L Banding performed on 20 metaphase cells	Expected karyotype	See Report
	Results: 46,XX Interpretation: This is a normal karyotype; no clonal abnormalities were detected at the stated band level of resolution.			
Post-Thaw Viable Cell Recovery	WiCell	Thaw using specified Thaw & Culture Recommendations	Recoverable attachment after passage	Pass
Identity by STR	WiCell	PowerPlex 16 HS System by Promega™	Defines STR profile of deposited cell line	See Report
Mycoplasma	WiCell	PCR	Amplification of mycoplasma specific DNA detected with negative result	Pass
Sterility	Steris	Native Product Direct Transfer using FTM and TSB (ST/07)	Negative for growth following 14 days of culture	Pass

Testing Reported by Provider

The Provider stated that some or all of the additional analyses listed below may have been performed for this cell line. For more information, publication and dbGaP links, where available, are provided on the cell line specific web page on the WiCell website.

- HumanCore Exome Kit
- Methylation
- Tra1-60 marker expression via flow cytometry
- Infinium® Expanded Multi-Ethnic Genotyping Array (MEGA^{EX})

Approval Date	WiCell Quality Assurance Approval
04-March-2025	3/4/2025 X HEB HEB WiCell Quality Assurance Signed by: HEBrunner

Date Reported: February 05, 2025

Cell Line: SCRP4301i-DB42912

Submitted Passage #: 14

Date of Sample: 1/30/2025


Specimen: Human IPSC

Results: 46,XX

Cell Line Sex: Female

Reason for Testing: LOT_RELEASE

Investigator: WiCell Stem Cell Bank, WiCell



Cell: 12

Slide: G03

Slide Type: Karyotype

Total Counted: 20

Total Analyzed: 8

Total Karyogrammed: 4

Band Resolution: 475 - 575

Interpretation:

This is a normal karyotype; no clonal abnormalities were detected at the stated band level of resolution.

Completed by: Erica Schutter, CG(ASCP)

Reviewed and Interpreted by: Justin Schleede, PhD, FACMG

For internal use only

Date: _____ **Sent By:** _____ **Sent To:** _____ **QC Review By:** _____

Limitations: This assay allows for microscopic visualization of numerical and structural chromosome abnormalities. The size of structural abnormality that can be detected is >3-10Mb, dependent upon the G-band resolution obtained from this specimen. For the purposes of this report, band level is defined as the number of G-bands per haploid genome. It is documented here as "band level", i.e., the range of bands determined from the four karyograms in this assay. Detection of heterogeneity of clonal cell populations in this specimen (i.e., mosaicism) is limited by the number of metaphase cells examined, documented here as "# of cells counted".

This assay was conducted solely for listed investigator/institution. The results of this assay are for research use only. Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.



Short Tandem Repeat

Requestor: WiCell Stem Cell Bank, WiCell

Sample Receipt Date: 30Jan25

STR Amplification Date: 06Feb25

Form SOP-89.01

Version 14.0

Sample Name	AI06e-SOX2YFP-WB68606 p68	SCR4301i-DB42912 p14	SCR4502i-DB42918 p14
WiCell CTR No. ¹	105671	105659	105658
FGA	Identifying information has been redacted to protect donor confidentiality. If more information is required, please contact info@wicell.org		
TPOX			
D8S1179			
vWA			
Amelogenin			
Penta_D			
CSF1PO			
D16S539			
D7S820			
D13S317			
D5S818			
Penta_E			
D18S51			
D21S11			
TH01			
D3S1358			
Allelic Polymorphisms	28	28	25
Matches ²	See Results		
Comments			

¹ CTR No.: Characterization Test Request Number; also known as a laboratory accessioning number.

² The STR profile of the sample(s) listed are a 100% match for the given sample unless otherwise specified.



Short Tandem Repeat

Requestor: WiCell Stem Cell Bank, WiCell

Sample Receipt Date: 30Jan25

STR Amplification Date: 06Feb25

Form SOP-89.01

Version 14.0

Assay Description: Short Tandem Repeat (STR) analysis is performed using the PowerPlex® 16 HS System by Promega™. Results are reported as 13 CODIS STR markers, Amelogenin for sex determination and two low-stutter, highly discriminating pentanucleotide STR markers.

Results: The genotypic profiles comprise a range of 25-28 allelic polymorphisms across the 15 STR loci analyzed. Sample 105671 is a 100% match to 105203, 104962, 104606, 104361, 103887, 103763, 101599, 99632, 98996, 98995 and additional profiles. Additional matches can be provided upon request.

Interpretation: The concentration of DNA required to achieve an acceptable STR genotype (signal/ noise) was equivalent to that required for the standard procedure (~1 ng/amplification reaction) from human genomic DNA. These results suggest that the cells submitted correspond to the cell lines as named and were not contaminated with any other human cells or a significant amount of mouse feeder layer cells.

Sensitivity: Sensitivity limits for detection of STR polymorphisms unique to either this or other human cell lines is ~2-4%.

2/13/2025	2/13/2025	2/17/2025
X Amber Kuhn	X Anna Lisa Larson	X Dawn Graham
Tech #1 Characterization Signed by: Kuhn, Amber	Tech #2 Characterization Signed by: Larson, Anna Lisa	QA Review Quality Assurance Signed by: Graham, Dawn

Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect. Raw data is available upon request.



Mycoplasma Assay Report

PCR-based assay performed by WiCell
WiCell Stem Cell Bank, WiCell
04Feb25

Form SOP-83.01
Version 6.0

Sample Name	Result	Interpretation
SCR4301i-DB42912 p15 (105739)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
SCR4206i-DB42864 p15 (105738)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
SCR4502i-DB42918 p15 (105708)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
WIBR3-S1-WB68685 p36 (105684)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
Positive (+) Control	Positive	
Negative (-) Control	Negative	

Assay Description

Sample is tested for presence of mycoplasma using EZ-PCR™ Mycoplasma Detection Kit (Sartorius).

2/4/2025	2/5/2025	2/5/2025
X Steph Dos Santos Tech #1 Characterization Signed by: Dos Santos, Stephany	X Kaylie Haddix Tech #2 Characterization Signed by: Haddix, Kaylie	X Dawn Graham QA Review Quality Assurance Signed by: Graham, Dawn

Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.

A gel image is available upon request.

Native Product Sterility Report



WiCell
504 S Rosa Road, Rm 101
Madison, WI 53719

SAMPLE #: 21041910
DATE RECEIVED: 29-Apr-21
TEST INITIATED: 03-May-21
TEST COMPLETED: 17-May-21

SAMPLE NAME / DESCRIPTION: WAI001-B-1-iETV2-WB67568
14746.001-WB67643
14747.001-WB67642
SCR4301i-DB42912
SCR4406i-DB42915
SCR4502i-DB42918
SCR4603i-DB42921
SCR4703i-DB42924
SCR4803i-DB42953
SCR4903i-DB42958

UNIQUE IDENTIFIER: N/A

TEST RESULTS:

# Tested	# Positives (Growth)	- Control
10	0	2 Negatives

TEST SUMMARY:

# Samples	Media Type	Volume (mL)	Incubation Temperature (° C)	Incubation Duration (Days)
10	TSB	40	20-25	14
10	FTG	40	30-35	14

REFERENCE: Processed according to LAB-003: Sterility Test Procedure

PD #: 000053

TEST METHODOLOGY: USP - Direct Transfer

COMMENTS: NA

REVIEWED BY Jonie Buckhand

DATE 17 May 2021

Specific test results may not be indicative of the characteristics of any other samples from the same lot or similar lots. This test report shall not be reproduced, except in full, without prior written approval. Liability is limited to the costs of the tests. Results applied to samples as received.