



Thaw and Culture Details

Cell Line Name	STAN251i-637C1
WiCell Lot Number	WB67455
Parent Material	STAN251i-637C1-DB44371
Provider	Stanford University – Laboratory of Dr. Thomas Quettermous
Banked By	WiCell
Thaw and Culture Recommendations	WiCell recommends thawing 1 vial into 3 wells of a 6 well plate.
Protocol	WiCell Feeder Independent mTeSR™ 1 Protocol
Culture Platform Prior to Freeze	Feeder Independent
	Medium: mTeSR™ 1
	Matrix: Matrigel®
Passage Number	p23 These cells were cultured for 22 passages prior to freeze and post colony selection. WiCell adds +1 to the passage number at freeze to best represent the overall passage number of the cells at thaw. Plated cells at thaw should be labeled passage 23.
Date Vialied	25-March-2020
Vial Label	STAN251i-637C1 p23 WB67455
Biosafety and Use Information	Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans.

Testing Performed by WiCell

Test Description	Test Provider	Test Method	Test Specification	Result
Karyotype by G-banding	WiCell	SOP-CH-003	Expected karyotype	See Report
Post-Thaw Viable Cell Recovery	WiCell	SOP-CH-305	≥ 15 Undifferentiated Colonies prior to passage, ≤ 30% Differentiation prior to passage, and recoverable attachment after passage	Pass
Identity by STR	UW Translational Research Initiatives in Pathology Laboratory	PowerPlex 16 HS System by Promega	Defines STR profile of deposited cell line	Pass
Sterility	Steris	ST/07	Negative	Pass
Mycoplasma	WiCell	SOP-CH-044	Negative	Pass

Testing Reported by Provider

The Provider stated that some or all of the additional analyses listed below may have been performed for this cell line. For more information, publication and dbGaP links, where available, are provided on the cell line specific web page on the WiCell website.

- RNA-Seq
- Whole Genome Sequencing



Approval Date	Quality Assurance Approval
02-July-2020	<div>7/17/2023</div> <div>X Ryen Smith</div> <div>JMG</div> <div>Quality Assurance</div> <div>Signed by: Smith, Ryen</div>

Date Reported: Monday, June 22, 2020

Cell Line: STAN251i-637C1-WB67455

Passage#: 23

Date of Sample: 6/15/2020


Specimen: Human IPSC

Results: 46,XX

Cell Line Sex: Female

Reason for Testing: LOT_RELEASE

Investigator: WiCell Stem Cell Bank, WiCell



Cell: 15

Slide: G03

Slide Type: Karyotype

Total Counted: 20

Total Analyzed: 8

Total Karyogrammed: 4

Band Resolution: 400 - 500

Interpretation:

This is a normal karyotype; no clonal abnormalities were detected at the stated band level of resolution.

Completed by: [REDACTED], CG(ASCP)

Reviewed and Interpreted by: [REDACTED] Ph.D.

Date: _____ **Sent By:** _____ **Sent To:** _____ **QC Review By:** _____

Limitations: This assay allows for microscopic visualization of numerical and structural chromosome abnormalities. The size of structural abnormality that can be detected is >3-10Mb, dependent upon the G-band resolution obtained from this specimen. For the purposes of this report, band level is defined as the number of G-bands per haploid genome. It is documented here as "band level", i.e., the range of bands determined from the four karyograms in this assay. Detection of heterogeneity of clonal cell populations in this specimen (i.e., mosaicism) is limited by the number of metaphase cells examined, documented here as "# of cells counted".

This assay was conducted solely for listed investigator/institution. The results of this assay are for research use only. Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.



Report Sent: 06/30/20

Short Tandem Repeat Analysis

Label on tube	STAN094i-081C2-WB67446 p.20 (81581)	STAN251i-637C1-WB67455 p.23 (81594)	CREM048i-BR3-1-WB67450 p.17 (81649)	STAN360i-465C2-WB67451 p.23 (81650)	GFAP-R416W-WB67486 p.18 (81714)	GFAP-R416R-WB67485 p.22 (81716)	MIN09i-33114.C-WB67492 p.17 (81717)
Label on Report	<p>Identifying information has been redacted to protect donor confidentiality. If more information is required, please, contact WiCell's Technical Support.</p>						
conc (ng/μL)							
A260/280							
Assay Date							
File Name							
FGA							
TPOX							
D8S1179							
vWA							
Amelogenin							
Penta_D							
CSF1PO							
D16S539							
D7S820							
D13S317							
D5S818							
Penta_E							
D18S51							
D21S11							
TH01							
D3S1358							
Allelic Polymorphisms							
Matches*							
Comments							



HISTOLOGY - IHC - MOLECULAR – IMAGING

Department of Pathology and Laboratory Medicine

TRIP Laboratory (Molecular)

<https://research.pathology.wisc.edu/trip-home/>

(608) 265-9168



Your Lab Partner

characterization@wicell.org

(608) 316-4145

Short Tandem Repeat Analysis

Label on tube	GFAP-R88R-WB67491 p.24 (81718)
Label on Report	Identifying information has been redacted to protect donor confidentiality. If more information is required, please, contact WiCell's Technical Support .
conc (ng/μL)	
A260/280	
Assay Date	
File Name	
FGA	
TPOX	
D8S1179	
vWA	
Amelogenin	
Penta_D	
CSF1PO	
D16S539	
D7S820	
D13S317	
D5S818	
Penta_E	
D18S51	
D21S11	
TH01	
D3S1358	
Allelic Polymorphisms	
Matches*	
Comments	

Short Tandem Repeat Analysis

Results: Based on the cells submitted by WiCell Characterization Department dated and received on 06/22/20, these samples define the STR profiles of the human cell lines as indicated by name. The genotypic profiles comprise a range of 23-27 allelic polymorphisms across the 15 STR loci analyzed.

Interpretation: The concentration of DNA required to achieve an acceptable STR genotype (signal/ noise) was equivalent to that required for the standard procedure (~1 ng/amplification reaction) from human genomic DNA. These results suggests that the cells submitted correspond to the cell lines as named and were not contaminated with any other human cells or a significant amount of mouse feeder layer cells.

Sensitivity: Sensitivity limits for detection of STR polymorphisms unique to either this or other human cell lines is ~2-5%.

Acknowledge TRIP in your publications, posters & presentations. For details, see:
<https://research.pathology.wisc.edu/acknowledging-trip/>

* **Note:** The STR profile of the following sample is an exact match for the given sample/samples.

X *RMB*

Digitally Signed on 06/30/20

_____, BA

TRIP Laboratory, Molecular

X *WMR*

Digitally Signed on 06/30/20

_____, PhD, Director / Co-Director

UWHC Molecular Diagnostics Laboratory / UWSMPH TRIP Laboratory

Testing was accomplished by analysis of human genetic polymorphisms at STR loci. This methodology has not yet been approved by the FDA and is for investigational use only.

Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at <https://www.wicell.org/media.acux/ca76d97c-862a-43f3-b02a-ab2d1e619100>. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.

Native Product Sterility Report



WiCell
504 S Rosa Road, Rm 101
Madison, WI 53719

SAMPLE #: 20051218
DATE RECEIVED: 21-May-20
TEST INITIATED: 28-May-20
TEST COMPLETED: 11-Jun-20

SAMPLE NAME / DESCRIPTION: WISCe011-A-40 WB67443
mCh-GnRH WB67447
GFAP-R88C WB67454
CREM024i-SS36-1 WB67440
H9-CAG-ChR2-EYFP WB67468
STAN094i-081C2 WB67446
MCW087i-U7112 WB67434
STAN360i-465C2 WB67451
CREM048i-BR3-1 WB67450
STAN251i-637C1 WB67455

UNIQUE IDENTIFIER: NA

TEST RESULTS:

# Tested	# Positives (Growth)	- Control
10	0	2 Negatives

TEST SUMMARY:

# Samples	Media Type	Volume (mL)	Incubation Temperature (° C)	Incubation Duration (Days)
10	TSB	40	20-25	14
10	FTG	40	30-35	14

REFERENCE: Processed according to LAB-003: Sterility Test Procedure

PD #: 000053

TEST METHODOLOGY: USP - Direct Transfer

COMMENTS: NA

REVIEWED BY

DATE 15 JUN 2020

Specific test results may not be indicative of the characteristics of any other samples from the same lot or similar lots. This test report shall not be reproduced, except in full, without prior written approval. Liability is limited to the costs of the tests. Results applied to samples as received.



Mycoplasma Assay Report

PCR-based assay performed by WiCell

WiCell

18Jun20

FORM SOP-CH-048.01

Version C Edition 01

Sample Name	Result	Comments/Suggestions
STAN094i-081C2-WB67446 (81578)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
CREM048i-BR3-1-WB67450 (81587)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
STAN251i-637C1-WB67455 (81588)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
STAN360i-465C2-WB67451 (81589)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
Positive (+) Control	Positive	
Negative (-) Control	Negative	

Reported by: [REDACTED] Cell Culture Specialist

Reviewed by: [REDACTED], Assistant Cell Culture Specialist

Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.

A gel image is available upon request.