

Cell Line	Steady-State Plasma Glucose Levels	Group	Clone
STAN303i-864C2 STAN304i-864C5	43	Insulin Sensitive	2 clones made from one donor.
STAN238i-553C5 STAN239i-553C8	48	Insulin Sensitive	2 clones made from one donor.
STAN269i-720C2 STAN270i-720C3	49	Insulin Sensitive	2 clones made from one donor.
STAN162i-345C2 STAN163i-345C5	50	Insulin Sensitive	2 clones made from one donor.
STAN131i-214C3 STAN132i-214C4	53	Insulin Sensitive	2 clones made from one donor.
STAN170i-364C3 STAN171i-364C4	54	Insulin Sensitive	2 clones made from one donor.
STAN226i-516C3 STAN227i-516C5	57	Insulin Sensitive	2 clones made from one donor.
STAN232i-541C2 STAN233i-541C5	60	Insulin Sensitive	2 clones made from one donor.
STAN276i-756C4 STAN277i-756C5	60	Insulin Sensitive	2 clones made from one donor.
STAN182i-410C3 STAN183i-410C4	61	Insulin Sensitive	2 clones made from one donor.
STAN197i-441C1 STAN198i-441C2	63	Insulin Sensitive	2 clones made from one donor.
STAN097i-107C3 STAN098i-107C4	66	Insulin Sensitive	2 clones made from one donor.
STAN261i-698C5	66	Insulin Sensitive	
STAN145i-286C1 STAN146i-286C2	67	Insulin Sensitive	2 clones made from one donor.
STAN187i-420C6	67	Insulin Sensitive	
STAN211i-475C1	69	Insulin Sensitive	
STAN222i-509C2 STAN223i-509C3	69	Insulin Sensitive	2 clones made from one donor.
STAN293i-835C1 STAN294i-835C5	70	Insulin Sensitive	2 clones made from one donor.
STAN166i-356C4 STAN167i-356C6	73	Insulin Sensitive	2 clones made from one donor.
STAN174i-371C3	73	Insulin Sensitive	
STAN315i-909C1 STAN316i-909C2	75	Insulin Sensitive	2 clones made from one donor.
STAN265i-707C1 STAN266i-707C4	76	Insulin Sensitive	2 clones made from one donor.
STAN301i-845AC2	76	Insulin Sensitive	2 clones made from one donor.

Cell Line	Steady-State Plasma Glucose Levels	Group	Clone
STAN302i-845AC4	76	Insulin Sensitive	2 clones made from one donor.
STAN091i-071C2	77	Insulin Sensitive	
STAN128i-211C2	80	Insulin Sensitive	
STAN295i-836C1	81	Insulin Sensitive	2 clones made from one donor.
STAN296i-836C2			
STAN284i-780C1	82	Insulin Sensitive	2 clones made from one donor.
STAN285i-780C2			
STAN331i-952C3	82	Insulin Sensitive	2 clones made from one donor.
STAN332i-952C5			
STAN333i-955C1	82	Insulin Sensitive	2 clones made from one donor.
STAN334i-955C2			
STAN338i-975C1	84	Insulin Sensitive	2 clones made from one donor.
STAN339i-975C3			
STAN154i-312C1	85	Insulin Sensitive	
STAN099i-108C2	88	Insulin Sensitive	2 clones made from one donor.
STAN100i-108C4			
STAN286i-799C2	89	Insulin Sensitive	
STAN317i-910C2	90	Insulin Sensitive	2 clones made from one donor.
STAN318i-910C3			
STAN350i-694C5	90	Insulin Sensitive	
STAN133i-215C1	91	Insulin Sensitive	2 clones made from one donor.
STAN134i-215C2			
STAN220i-504C2	92	Insulin Sensitive	
STAN206i-459C1	93	Insulin Sensitive	2 clones made from one donor.
STAN207i-459C2			
STAN148i-291C3	94	Insulin Sensitive	
STAN135i-230C1	95	Insulin Sensitive	2 clones made from one donor.
STAN136i-230C2			
STAN114i-135C2	98	Insulin Sensitive	
STAN164i-352C1	100	Insulin Sensitive	2 clones made from one donor.
STAN165i-352C5			
STAN341i-994C1	102	Insulin Sensitive	2 clones made from one donor.
STAN342i-994C3			
STAN343i-998C1	102	Insulin Sensitive	
STAN345i-390C3	102	Insulin Sensitive	
STAN119i-192C1	103	Insulin Sensitive	2 clones made from one donor.
STAN120i-192C2			
STAN313i-907C1	103	Insulin Sensitive	2 clones made from one donor.
STAN314i-907C4			

Cell Line	Steady-State Plasma Glucose Levels	Group	Clone
STAN327i-940C2	103	Insulin Sensitive	
STAN202i-445C1 STAN203i-445C2	104	Insulin Sensitive	2 clones made from one donor.
STAN248i-617C1 STAN249i-617C2	105	Insulin Sensitive	2 clones made from one donor.
STAN255i-649C1 STAN256i-649C2	106	Insulin Sensitive	2 clones made from one donor.
STAN287i-810C1 STAN288i-810C3	106	Insulin Sensitive	2 clones made from one donor.
STAN084i-041C1 STAN085i-041C3	108	Intermediate	2 clones made from one donor.
STAN340i-981C1	112	Intermediate	
STAN155i-313C2	114	Intermediate	
STAN156i-334C1 STAN157i-334C2	115	Intermediate	2 clones made from one donor.
STAN247i-614C2	117	Intermediate	
STAN212i-488C2 STAN213i-488C3	120	Intermediate	2 clones made from one donor.
STAN104i-120C3	129	Intermediate	
STAN191i-429C2 STAN192i-429C3	130	Intermediate	2 clones made from one donor.
STAN257i-685C3 STAN258i-685C4	132	Intermediate	2 clones made from one donor.
STAN179i-400C2	133	Intermediate	
STAN095i-102C4 STAN096i-102C6	136	Intermediate	2 clones made from one donor.
STAN311i-906C1 STAN312i-906C3	136	Intermediate	2 clones made from one donor.
STAN210i-468C1	138	Intermediate	
STAN280i-770C1 STAN281i-770C2	138	Intermediate	2 clones made from one donor.
STAN129i-212C2 STAN130i-212C4	141	Intermediate	2 clones made from one donor.
STAN080i-023C2 STAN081i-023C3	142	Intermediate	2 clones made from one donor.
STAN218i-501C5 STAN219i-501C6	142	Intermediate	2 clones made from one donor.
STAN086i-054C1 STAN087i-054C2	145	Intermediate	2 clones made from one donor.
STAN149i-298C1	145	Intermediate	

Cell Line	Steady-State Plasma Glucose Levels	Group	Clone
STAN150i-337C1	145	Intermediate	
STAN204i-448C1 STAN205i-448C2	145	Intermediate	2 clones made from one donor.
STAN144i-285C3	147	Intermediate	
STAN082i-024C1 STAN083i-024C4	153	Intermediate	2 clones made from one donor.
STAN299i-843C5 STAN300i-843C6	153	Intermediate	2 clones made from one donor.
STAN193i-432C1 STAN194i-432C3	157	Intermediate	2 clones made from one donor.
STAN245i-601C4 STAN246i-601C5	161	Intermediate	2 clones made from one donor.
STAN151i-303C3	163	Intermediate	
STAN208i-460C2 STAN209i-460C6	163	Intermediate	2 clones made from one donor.
STAN175i-373C4 STAN176i-373C6	165	Intermediate	2 clones made from one donor.
STAN251-637C1 STAN252i-637C2	168	Intermediate	2 clones made from one donor.
STAN093i-081C1 STAN094i-081C2	171	Intermediate	2 clones made from one donor.
STAN215i-490C3	171	Intermediate	
STAN240i-558C3 STAN241i-558C4	171	Intermediate	2 clones made from one donor.
STAN321i-924C1 STAN322i-924C2	172	Intermediate	2 clones made from one donor.
STAN224i-514C3 STAN225i-514C4	173	Intermediate	2 clones made from one donor.
STAN282i-771C3 STAN283i-771C4	173	Intermediate	2 clones made from one donor.
STAN108i-123C2 STAN109i-123C3	174	Intermediate	2 clones made from one donor.
STAN160i-341C4 STAN161i-341C5	175	Intermediate	2 clones made from one donor.
STAN250i-622C2	175	Intermediate	
STAN278i-767C2 STAN279i-767C3	178	Intermediate	2 clones made from one donor.
STAN078i-018C2 STAN079i-018C3	180	Intermediate	2 clones made from one donor.
STAN125i-202C2 STAN126i-202C3	181	Intermediate	2 clones made from one donor.

Cell Line	Steady-State Plasma Glucose Levels	Group	Clone
STAN291i-827C1 STAN292i-827C2	181	Intermediate	2 clones made from one donor.
STAN111i-125C4	186	Intermediate	
STAN141i-250C4	187	Intermediate	
STAN230i-533C1 STAN231i-533C2	193	Insulin Resistant	2 clones made from one donor.
STAN088i-060C1	196	Insulin Resistant	
STAN076i-004C1 STAN077i-004C2	197	Insulin Resistant	2 clones made from one donor.
STAN152i-307C2 STAN153i-307C3	199	Insulin Resistant	2 clones made from one donor.
STAN216i-496C1 STAN217i-496C2	199	Insulin Resistant	2 clones made from one donor.
STAN195i-434C3	200	Insulin Resistant	
STAN297i-839C1 STAN298i-839C5	200	Insulin Resistant	2 clones made from one donor.
STAN323i-928C2 STAN324i-928C6	201	Insulin Resistant	2 clones made from one donor.
STAN200i-444C1 STAN201i-444C2	202	Insulin Resistant	2 clones made from one donor.
STAN234i-546C1 STAN235i-546C2	202	Insulin Resistant	2 clones made from one donor.
STAN307i-871C2 STAN308i-871C3	202	Insulin Resistant	2 clones made from one donor.
STAN319i-922C1	207	Insulin Resistant	
STAN329i-949C3 STAN330i-949C4	208	Insulin Resistant	2 clones made from one donor.
STAN140i-243C1	209	Insulin Resistant	
STAN168i-358C1 STAN169i-358C4	212	Insulin Resistant	2 clones made from one donor.
STAN328i-946C2	212	Insulin Resistant	
STAN228i-523C1 STAN229i-523C2	215	Insulin Resistant	2 clones made from one donor.
STAN173i-368C2	217	Insulin Resistant	
STAN263i-703C1 STAN264i-703C2	231	Insulin Resistant	2 clones made from one donor.
STAN106i-121C1 STAN107i-121C2	233	Insulin Resistant	2 clones made from one donor.
STAN253i-647C1	235	Insulin Resistant	2 clones made from one donor.

Cell Line	Steady-State Plasma Glucose Levels	Group	Clone
STAN254i-647C3	235	Insulin Resistant	2 clones made from one donor.
STAN305i-865C2 STAN306i-865C3	235	Insulin Resistant	2 clones made from one donor.

Cell Line	Steady-State Plasma Glucose Levels	Group	Clone
STAN190i-424C2	236	Insulin Resistant	
STAN158i-336C2 STAN159i-336C3	243	Insulin Resistant	2 clones made from one donor.
STAN273i-729C1 STAN274i-729C2	244	Insulin Resistant	2 clones made from one donor.
STAN137i-231C3	245	Insulin Resistant	
STAN275i-732C1	248	Insulin Resistant	
STAN271i-726C3 STAN272i-726C4	251	Insulin Resistant	2 clones made from one donor.
STAN180i-408C2 STAN181i-408C4	253	Insulin Resistant	2 clones made from one donor.
STAN121i-193C2 STAN122i-193C1	255	Insulin Resistant	2 clones made from one donor.
STAN242i-573C2	257	Insulin Resistant	
STAN336i-963C1 STAN337i-963C3	258	Insulin Resistant	2 clones made from one donor.
STAN243i-589C3 STAN244i-589C4	260	Insulin Resistant	2 clones made from one donor.
STAN115i-148C2 STAN116i-148C1	266	Insulin Resistant	2 clones made from one donor.
STAN089i-070C2 STAN090i-070C4	267	Insulin Resistant	2 clones made from one donor.
STAN138i-239C3 STAN139i-239C5	272	Insulin Resistant	2 clones made from one donor.
STAN260i-688C3	279	Insulin Resistant	
STAN117i-182C4 STAN118i-182C7	280	Insulin Resistant	2 clones made from one donor.
STAN236i-551C1 STAN237i-551C2	282	Insulin Resistant	2 clones made from one donor.
STAN142i-255C1	323	Insulin Resistant	
STAN346i-59C2	Not provided	Not provided	
STAN347i-66C3	Not provided	Not provided	
STAN348i-851C1	Not provided	Not provided	
STAN349i-762C3	Not provided	Not provided	
STAN351i-631C1	Not provided	Not provided	
STAN357i-298C2 STAN358i-298C3	Not provided	Not provided	2 clones made from one donor.
STAN359i-442C11	Not provided	Not provided	
STAN360i-465C2	Not provided	Not provided	

Cell Line	Steady-State Plasma Glucose Levels	Group	Clone
STAN361i-531C6	Not provided	Not provided	
STAN362i-178C1	Not provided	Not provided	
STAN363i-845C5	Not provided	Not provided	
STAN364i-055C2	Not provided	Not provided	
STAN365i-664C1	Not provided	Not provided	
STAN366i-282C2	Not provided	Not provided	
STAN367i-894C1	Not provided	Not provided	2 clones made from one donor.
STAN368i-894C2			
STAN369i-088C3	Not provided	Not provided	
STAN370i-190C3	Not provided	Not provided	
STAN371i-868C5	Not provided	Not provided	
STAN372i-411C8	Not provided	Not provided	
STAN373i-412C4	Not provided	Not provided	
STAN374i-565C4	Not provided	Not provided	
STAN375i-961C1	Not provided	Not provided	
STAN376i-518C2	Not provided	Not provided	
STAN377i-572C5	Not provided	Not provided	
STAN378i-886C4	Not provided	Not provided	
STAN379i-783C6	Not provided	Not provided	
STAN380i-092C2	Not provided	Not provided	
STAN381i-652C1	Not provided	Not provided	