



## Thaw and Culture Details


Cell Line Name	CREM041i-SA208-2
WiCell Lot Number	DB48099
Provider	Boston University – Laboratory of Dr. Martin Steinberg
Banked By	Boston University - Laboratory of Dr. Gustavo Mostoslavsky
Thaw and Culture Recommendations	The Provider recommends thawing 1 vial into 2 wells of a 6 well plate.
Culture Platform	Feeder Dependent
	Medium: hESC Medium (KOSR)
	Matrix: MEF
Protocol	WiCell Feeder Dependent Protocol
Passage Number	p6 These cells were cultured for 6 passages after colony picking prior to freeze. Add +1 to the passage number to best represent the overall passage number of the cells at thaw.
Date Viald	31-May-2016
Vial Label	SA208-2p6 hiPSC/KSR 5/31/16 SMP
Biosafety and Use Information	Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans.

## Testing Reported by Provider

The Provider stated that some or all of the additional analyses listed below may have been performed for this cell line. For more information, publication and dbGaP links, where available, are provided on the cell line specific web page on the WiCell website.

- Digital Genome Sequencing
- Infinium® Expanded Multi-Ethnic Genotyping Array (MEGA<sup>EX</sup>)

Please note: Prior to shipment of these cells, WiCell will perform the following characterization assays: post-thaw viable recovery, identity by STR, sterility, mycoplasma, and karyotype.

Approval Date	Quality Assurance Approval
05-December-2016	<div style="text-align: right;">12/5/2016</div> <div style="text-align: center;"> X AMK AMK Quality Assurance Signed by Klade, Anjelica</div>