



Thaw and Culture Details

Cell Line Name	CREM061i-BT1-1
WiCell Lot Number	DB66780
Provider	Boston University – Laboratory of Dr. Martin Steinberg
Banked By	Boston University - Laboratory of Dr. Gustavo Mostoslavsky
Thaw and Culture Recommendations	The Provider recommends thawing 1 vial into 2 wells of a 6 well plate.
Culture Platform	Feeder Independent
	Medium: mTeSR™1
	Matrix: Matrigel®
Protocol	WiCell Feeder Independent mTeSR™1 Protocol
Passage Number	p25 These cells were cultured for 25 passages prior to freeze and post colony picking. Therefore, plated cells at thaw should be labeled passage 26.
Date Vialled	12-March-2015
Vial Label	B-THAL1-1p25 hiPSC/mTeSR 3/11/15 SMP
Biosafety and Use Information	Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans.

Testing Reported by Provider

The Provider stated that some or all of the additional analyses listed below may have been performed for this cell line. For more information, publication and dbGaP links, where available, are provided on the cell line specific web page on the WiCell website.

- Digital Genome Sequencing
- Infinium® Expanded Multi-Ethnic Genotyping Array (MEGA^{EX})

Please note: Prior to shipment of these cells, WiCell will perform the following characterization assays: post-thaw viable recovery, identity by STR, sterility, mycoplasma, and karyotype.

Approval Date	Quality Assurance Approval
25-June-2018	<div style="text-align: right; font-size: small;">6/25/2018</div> <div style="text-align: center;">X HEB HEB Quality Assurance Signed by Bruner, Haley</div>