



Thaw and Culture Details

| | |
|----------------------------------|---|
| Cell Line Name | JHU095i |
| WiCell Lot Number | DB41261 |
| Provider | Johns Hopkins University – Laboratory of Dr. Lewis Becker |
| Banked By | Johns Hopkins University – Laboratory of Dr. Lewis Becker |
| Thaw and Culture Recommendations | The Provider recommends thawing 1 vial into 2 wells of a 6 well plate. |
| Culture Platform | Feeder Independent |
| | Medium: E8 |
| | Matrix: Vitronectin |
| Protocol | WiCell Feeder Independent E8 Medium Protocol |
| Passage Number | p5 These cells were cultured for 5 passages post reprogramming prior to freeze. Add +1 to the passage number to best represent the overall passage number of the cells at thaw. |
| Date Viald | 20-May-2016 |
| Vial Label | P95 P5 0.7X10 ⁶ 5/20/16 |
| Biosafety and Use Information | This cell line is of human origin. Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans. |

Testing Reported by Provider

The Provider stated that some or all of the additional analyses listed below may have been performed for this cell line. For more information, publication and dbGaP links, where available, are provided on the cell line specific web page on the WiCell website.

- Embryoid bodies
- Infinium® Expanded Multi-Ethnic Genotyping Array (MEGA^{EX})

Please note: Prior to shipment of these cells, WiCell will perform the following characterization assays: post-thaw viable recovery, identity by STR, sterility, mycoplasma, and karyotype.

| Approval Date | Quality Assurance Approval |
|----------------|----------------------------|
| 26-August-2016 | |