



Thaw and Culture Details

Cell Line Name	MCW004i-40002545
WiCell Lot Number	WB66488
Provider	Medical College of Wisconsin – Laboratory of Dr. Ulrich Broeckel
Banked By	WiCell
Thaw and Culture Recommendations	WiCell recommends thawing 1 vial into 4 wells of a 6 well plate.
Culture Platform	Feeder Independent
	Medium: mTeSR™1
	Matrix: Matrigel®
Protocol	WiCell Feeder Independent mTeSR1 Protocol
Passage Number	p13 These cells were cultured for 12 passages prior to freeze and post colony picking. WiCell adds +1 to the passage number at freeze to best represent what the overall passage number of the cells at thaw. Plated cells at thaw should be labeled passage 13.
Date Vialied	04-August-2017
Vial Label	MCW004i-40002545 p13 WB66488
Biosafety and Use Information	Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans.

Testing Performed by WiCell

Test Description	Test Provider	Test Method	Test Specification	Result
Karyotype by G-banding	WiCell	SOP-CH-003	Expected karyotype	See Report
	Results: 47,XX+8[3]/46,XX[17] Interpretation: This is an abnormal karyotype; results contain a clonal chromosomal aberration recurrently acquired in cultures of this cell type. No other clonal abnormalities were detected at the stated band level of resolution.			
Post-Thaw Viable Cell Recovery	WiCell	SOP-CH-305	≥ 15 Undifferentiated Colonies prior to passage, ≤ 30% Differentiation prior to passage, and recoverable attachment after passage	Pass
Identity by STR	UW Translational Research Initiatives in Pathology Laboratory	PowerPlex 16 HS System by Promega	Defines STR profile of deposited cell line	Pass
Sterility	Steris	ST/07	Negative	Pass
Mycoplasma	WiCell	SOP-QU-004	Negative	Pass



Testing Reported by Provider

The Provider stated that some or all of the additional analyses listed below may have been performed for this cell line. For more information, publication and dbGaP links, where available, are provided on the cell line specific web page on the WiCell website.

- Tra1-60 marker expression
- mRNA expression by qPCR
- Infinium® Expanded Multi-Ethnic Genotyping Array (MEGA^{EX})

Approval Date	Quality Assurance Approval
14-May-2018	<p style="text-align: right;">5/21/2020</p> <p>X HEB _____ HEB Quality Assurance Signed by: Bruner, Haley</p>

Date Reported: Monday, June 04, 2018
Cell Line: MCW004i-40002545-WB66488
13751

Cell Line Sex: Female
Reason for Testing: lot release testing


Passage#: 13

Date of Sample: 5/25/2018

Investigator: [REDACTED] WiCell

Specimen: Human IPS

Results: 47,XX+8[3]/46,XX[17]



Cell: 28
Slide: G03
Slide Type: Karyotype

Total Counted: 20
Total Analyzed: 9
Total Karyogrammed: 5
Band Resolution: 375 - 500

Interpretation:

This is an abnormal karyotype; results contain a clonal chromosomal aberration recurrently acquired in cultures of this cell type. No other clonal abnormalities were detected at the stated band level of resolution.

Completed by: [REDACTED]

Reviewed and Interpreted by: [REDACTED], PhD, FACMG

A signed copy of this report is available upon request.

Date: _____ **Sent By:** _____ **Sent To:** _____ **QC Review By:** _____

Limitations: This assay allows for microscopic visualization of numerical and structural chromosome abnormalities. The size of structural abnormality that can be detected is >3-10Mb, dependent upon the G-band resolution obtained from this specimen. For the purposes of this report, band level is defined as the number of G-bands per haploid genome. It is documented here as "band level", i.e., the range of bands determined from the four karyograms in this assay. Detection of heterogeneity of clonal cell populations in this specimen (i.e., mosaicism) is limited by the number of metaphase cells examined, documented here as "# of cells counted".

This assay was conducted solely for listed investigator/institution. The results of this assay are for research use only. Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.

Sample Report:

13751-STR

Sample Name on Tube: 13751-STR91.6 ng/ μ L, (A260/280=1.89)**Sample Type:** Cells**Cell Count:** ~2 million cells**Requestor:**

WiCell Research Institute

Quality Department

Sample Date: N/A**Receive Date:** 05/29/18**Assay Date:** 05/29/18**File Name:** STR 180530c wmr**Report Date:** 06/04/18

STR Locus	STR Genotype Repeat #	STR Genotype
FGA	16-18,18.2,19,19.2,20,20.2,21,21.2,22, 22.2, 23, 23.2, 24, 24.2, 25, 25.2, 26-30, 31.2, 43.2, 44.2,45.2, 46.2	Identifying information has been redacted to protect donor confidentiality. If more information is required, please, contact WiCell's Technical Support .
TPOX	6-13	
D8S1179	7-18	
vWA	10-22	
Amelogenin	X,Y	
Penta D	2.2, 3.2, 5, 7-17	
CSF1PO	6-15	
D16S539	5, 8-15	
D7S820	6-14	
D13S317	7-15	
D5S818	7-16	
Penta E	5-24	
D18S51	8-10, 10.2, 11-13, 13.2, 14-27	
D21S11	24,24.2,25,25.2,26-28,28.2,29,29.2, 30, 30.2,31, 31.2,32,32.2,33,33.2, 34,34.2,35,35.2,36-38	
TH01	4-9,9.3,10-11,13.3	
D3S1358	12-20	

Results: Based on the 13751-STR cells submitted by WiCell QA dated and received on 05/29/18, this sample (Label on Tube: 13751-STR) defines the STR profile of the human stem cell line MCW004i-40002545 comprising 28 allelic polymorphisms across the 15 STR loci analyzed.

Interpretation: No STR polymorphisms other than those corresponding to the human MCW004i-40002545 stem cell line were detected and the concentration of DNA required to achieve an acceptable STR genotype (signal/ noise) was equivalent to that required for the standard procedure (~1 ng/amplification reaction) from human genomic DNA. This result suggests that the 13751-STR sample submitted corresponds to the MCW004i-40002545 stem cell line and was not contaminated with any other human stem cells or a significant amount of mouse feeder layer cells.

Sensitivity: Sensitivity limits for detection of STR polymorphisms unique to either this or other human stem cell lines is ~2-5%.



Digitally Signed on 06/05/18

_____, BA
TRIP Laboratory, Molecular

Digitally Signed on 06/05/18

_____, PhD, Director / Co-Director
UWHC Molecular Diagnostics Laboratory / UWSMPH TRIP Laboratory

Native Product Sterility Report



WiCell
504 S Rosa Rd, Rm 101
Madison, WI 53719

SAMPLE #: 18030537
DATE RECEIVED: 08-Mar-18
TEST INITIATED: 13-Mar-18
TEST COMPLETED: 27-Mar-18

SAMPLE NAME / DESCRIPTION: STAN129i-212C2 DB35772 13516
STAN130i-212C4 DB35777 13517
MCW002i-40001265 WB66495 13518
MCW004i-40002545 WB66488 13519
MCW006i-40000930 WB66499 13520
MCW008i-40000992 WB66496 13521
MCW010i-40000756 WB66487 13522
MCW011i-40000664 WB66486 13523
MCW015i-A2196 WB66497 13524
MCW016i-A2159 WB66510 13525
MCW021i-50001743 WB66448 13526
MCW025i-A2566 WB66504 13527
MCW034i-A2780 WB66502 13528
MCW036i-A3170 WB66501 13529
MCW037i-50000777 WB66459 13530
MCW041i-U2104 WB66494 13531
MCW048i-40001845 WB66460 13532
MCW050i-40000626 WB66467 13533
MCW067i-40001036 WB66478 13534
MCW068i-40002385 WB66452 13535

UNIQUE IDENTIFIER: NA
PRODUCT REGISTRATION: Other: Human iPS cells

TEST RESULTS:

# Tested	# Positives (Growth)	- Control
20	0	2 Negatives

TEST SUMMARY:

# Samples	Media Type	Volume (mL)	Incubation Temperature (° C)	Incubation Duration (Days)
20	TSB	40	20-25	14
20	FTG	40	30-35	14

REFERENCE: Processed according to LAB-003: Sterility Test Procedure

METHOD VALIDATION / PD #: 000053

Native Product Sterility Report



TEST METHODOLOGY: USP - Direct Transfer

COMMENTS: Sample #18030537
Report as per packing slip.

REVIEWED BY *Olesand*

DATE 28MAR18

Specific test results may not be indicative of the characteristics of any other samples from the same lot or similar lots. This test report shall not be reproduced, except in full, without prior written approval. Liability is limited to the costs of the tests.



Mycoplasma Detection Assay Report

Testing Performed by WiCell

Lot Release Testing

May 24, 2018

FORM SOP-QU-004.01

Version G Edition 02

Reported by: AP

Reviewed by: DF

BD Monolight 180

#	Sample Name	Reading A		A	Reading B		B	Ratio B/A	Result	Comments/Suggestions
		RLU1	RLU2	Ave	RLU1	RLU2	Ave			
1	MCW004i-40002545-WB66488 13751	325	294	309.5	124	121	122.5	0.40	Negative	
2	Positive (+) Control	605	627	616	33968	34195	34082	55.33	Positive	
3	Negative (-) Control	1030	1010	1020	105	117	111	0.11	Negative	

