




Thaw and Culture Details

Cell Line Name	H1-FMR1-KO
WiCell Lot Number	WB67517
Provider	University of Wisconsin - Dr. Anita Bhattacharyya and Dr. Xinyu Zhao
Banked By	WiCell
Thaw and Culture Recommendations	WiCell recommends thawing 1 vial into 3 wells of a 6 well plate.
Culture Platform	Feeder Independent
	Medium: mTeSR™1
	Matrix: Matrigel®
Protocol	WiCell Feeder Independent mTeSR™1 Protocol
Passage Number	p40 These cells were cultured for 39 passages prior to freeze. Cells were modified from WA01 at passage 30. WiCell adds +1 to the passage number at freeze to best represent the overall passage number of the cells at thaw. Plated cells at thaw should be labeled passage 40.
Date Vialied	18-July-2020
Vial Label	H1-FMR1-KO p40 WB67517
Biosafety and Use Information	Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans.

Testing Performed by WiCell

Test Description	Test Provider	Test Method	Test Specification	Result
Karyotype by G-banding	WiCell	SOP-CH-003	Expected karyotype	See Report
Post-Thaw Viable Cell Recovery	WiCell	SOP-CH-305	≥ 15 Undifferentiated Colonies prior to passage, ≤ 30% Differentiation prior to passage, and recoverable attachment after passage	Pass
Identity by STR	UW Translational Research Initiatives in Pathology Laboratory	PowerPlex 16 HS System by Promega	Defines STR profile of deposited cell line	Pass
Sterility	Steris	ST/07	Negative	Pass
Mycoplasma	WiCell	SOP-CH-044	Negative	Pass

Approval Date	Quality Assurance Approval
27-August-2020	<div style="text-align: right;">8/27/2020</div>  <small>X JKG JKG Quality Assurance Signed by Gay, Jenna</small>

Date Reported: Friday, August 7, 2020

Cell Line: H1-FMR1-KO-WB67517

Submitted Passage #: 40

Date of Sample: 8/3/2020


Specimen: Human Modified ESC

Results: 46,XY

Cell Line Sex: Male

Reason for Testing: LOT_RELEASE

Investigator: WiCell Stem Cell Bank, WiCell



Cell: 1

Slide: G01

Slide Type: Karyotype

Total Counted: 20

Total Analyzed: 8

Total Karyogrammed: 4

Band Resolution: 425 - 475

Interpretation:

This is a normal karyotype; no clonal abnormalities were detected at the stated band level of resolution.

Completed by: [REDACTED], CG(ASCP)

Reviewed and Interpreted by: [REDACTED], Ph.D.

Date: _____ **Sent By:** _____ **Sent To:** _____ **QC Review By:** _____

Limitations: This assay allows for microscopic visualization of numerical and structural chromosome abnormalities. The size of structural abnormality that can be detected is >3-10Mb, dependent upon the G-band resolution obtained from this specimen. For the purposes of this report, band level is defined as the number of G-bands per haploid genome. It is documented here as "band level", i.e., the range of bands determined from the four karyograms in this assay. Detection of heterogeneity of clonal cell populations in this specimen (i.e., mosaicism) is limited by the number of metaphase cells examined, documented here as "# of cells counted".

This assay was conducted solely for listed investigator/institution. The results of this assay are for research use only. Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.

Short Tandem Repeat Analysis

Receive Date: 08/05/20

Report Sent: 08/12/20

Requestor: WiCell Characterization

Label on tube	82127	82128	82129	82154	82181	82204	82205
Label on Report	JHU027i-DB40972 p.7 (82127)	H1-FMR1-FLAG- WB67514 p.36 (82128)	STAN130i-212C4- WB67515 p.16 (82129)	CREM005i-SS2-1GAG- DB66769 p.42 (82154)	STAN120i-192C2- WB67516 p.18 (82181)	H1-FMR1-KO-WB67517 p.40 (82204)	STAN099i-108C2- WB67518 p.15 (82205)
conc (ng/μL)	Identifying information has been redacted to protect donor confidentiality. If more information is required, please, contact WiCell's Technical Support.					129.5	Identifying information has been redacted to protect donor confidentiality. If more information is required, please, contact WiCell's Technical Support.
A260/280						1.82	
Assay Date						8/6/2020	
File Name						STR 200807 wmr	
FGA						20,24	
TPOX						8,11	
D8S1179						12,13	
vWA						15,17	
Amelogenin						X,Y	
Penta_D						10,13	
CSF1PO						12,13	
D16S539						9,13	
D7S820						8,12	
D13S317						8,11	
D5S818						9,11	
Penta_E						10,12	
D18S51						17,18	
D21S11	28,32.2						
TH01	9.3,9.3						
D3S1358	15,15						
Allelic Polymorphisms	27	28	28	27	25	28	26
Matches*			80512, 70862		77321		79403
Comments							

Short Tandem Repeat Analysis

Label on tube	82206
Label on Report	STAN378i-886C4-WB67520 p.27 (82206)
conc (ng/μL)	Identifying information has been redacted to protect donor confidentiality. If more information is required, please, contact WiCell's Technical Support .
A260/280	
Assay Date	
File Name	
FGA	
TPOX	
D8S1179	
vWA	
Amelogenin	
Penta_D	
CSF1PO	
D16S539	
D7S820	
D13S317	
D5S818	
Penta_E	
D18S51	
D21S11	
TH01	
D3S1358	
Allelic Polymorphisms	28
Matches*	77678
Comments	



HISTOLOGY - IHC - MOLECULAR – IMAGING
 Department of Pathology and Laboratory Medicine
 TRIP Laboratory (Molecular)
<https://research.pathology.wisc.edu/trip-home/>
 (608) 265-9168



Your Lab Partner
 characterization@wicell.org
 (608) 316-4145

Short Tandem Repeat Analysis

Results: Based on the DNA submitted by WiCell Characterization Department dated and received on 08/05/20, these samples define the STR profiles of the human cell lines as indicated by name. The genotypic profiles comprise a range of 25-28 allelic polymorphisms across the 15 STR loci analyzed.

Interpretation: The concentration of DNA required to achieve an acceptable STR genotype (signal/ noise) was equivalent to that required for the standard procedure (~1 ng/amplification reaction) from human genomic DNA. These results suggests that the cells submitted correspond to the cell lines as named and were not contaminated with any other human cells or a significant amount of mouse feeder layer cells.

Sensitivity: Sensitivity limits for detection of STR polymorphisms unique to either this or other human cell lines is ~2-5%.

**Acknowledge TRIP in your publications, posters & presentations. For details, see:
<https://research.pathology.wisc.edu/acknowledging-trip/>**

* **Note:** The STR profile of the following sample is an exact match for the given sample/samples.

X RMB Digitally Signed on **08/12/20**
 _____, BA
 TRIP Laboratory, Molecular

X WMR Digitally Signed on **08/12/20**
 _____, PhD, Director / Co-Director
 UWHC Molecular Diagnostics Laboratory / UWSMPH TRIP Laboratory

Testing was accomplished by analysis of human genetic polymorphisms at STR loci. This methodology has not yet been approved by the FDA and is for investigational use only.
 Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at <https://www.wicell.org/media.acux/ca76d97c-862a-43f3-b02a-ab2d1e619100>. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.

Native Product Sterility Report



WiCell
504 S Rosa Road, Rm 101
Madison, WI 53719

SAMPLE #: 20071581
DATE RECEIVED: 30-Jul-20
TEST INITIATED: 31-Jul-20
TEST COMPLETED: 14-Aug-20

SAMPLE NAME / DESCRIPTION: MIN12i-33362.C WB67499
WISCe011-A-39 WB67500
STAN120i-192C2 WB67516
STAN130i-212C4 WB67515
STAN099i-108C2 WB67518
H1-FMR1-FLAG WB67514
H1-FMR1-KO WB67517
STAN378i-886C4 WB67520
STAN206i-459C1 WB67519
JHU027i DB40972

UNIQUE IDENTIFIER: NA

TEST RESULTS:

# Tested	# Positives (Growth)	- Control
10	0	2 Negatives

TEST SUMMARY:

# Samples	Media Type	Volume (mL)	Incubation Temperature (° C)	Incubation Duration (Days)
10	TSB	40	20-25	14
10	FTG	40	30-35	14

REFERENCE: Processed according to LAB-003: Sterility Test Procedure

PD #: 000053

TEST METHODOLOGY: USP - Direct Transfer

COMMENTS: NA

REVIEWED BY

DATE 24 Aug 2020

Native Product Sterility Report



Specific test results may not be indicative of the characteristics of any other samples from the same lot or similar lots. This test report shall not be reproduced, except in full, without prior written approval. Liability is limited to the costs of the tests. Results applied to samples as received.



Mycoplasma Assay Report

PCR-based assay performed by WiCell

WiCell

04Aug20

FORM SOP-CH-048.01

Version C Edition 01

Sample Name	Result	Comments/Suggestions
STAN206i-459C1-WB67519 (82200)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
STAN378i-886C4-WB67520 (82201)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
STAN099i-108C2-WB67518 (82202)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
H1-FMR1-KO-WB67517 (82203)	Negative	Band was not seen at 270bp, indicating the absence of mycoplasma.
Positive (+) Control	Positive	
Negative (-) Control	Negative	

Reported by: [REDACTED], Research Specialist

Reviewed by: [REDACTED], Assistant Cell Culture Specialist

Unless otherwise mutually agreed in writing, the services provided to you hereunder by WiCell Research Institute, Inc. ("WiCell") are governed solely by WiCell's Terms and Conditions of Service, found at www.wicell.org/privacyandterms. Any terms you may attach to a purchase order or other document that are inconsistent, add to, or conflict with WiCell's Terms and Conditions of Service are null and void and of no legal force or effect.

A gel image is available upon request.